


NURSE COURIER® 3
MARQUEE
Software



Critical Information at a Glance

Telecor's Marquee Software provides a link between the Telecor Nurse Courier System and digital message displays. This versatile application allows active calls in the nurse call system to be shown on remote displays. Conveniently located throughout your facility, these displays give your organization wider access to view and monitor call activity on the system.

high, medium and low priorities



Calls are displayed in priority sequence by call type using the scrolling effect on the display. Call types can be represented in varying colors, making priority calls easily distinguishable from standard calls. When there are no active calls, a user-defined description for the display (such as the location) scrolls across the screen, in addition to the current date and time.

Displays also provide an excellent means of alerting a specific location to particular types of calls. For example, there may be staff locations where personnel need to be notified of any Code Blue calls that occur, system wide. The Marquee software can filter calls at any selected display to only show calls with a specified status.

Installation is easy and inexpensive. The displays can reside on the facility's Local Area Network (LAN) so they can be installed at any location within the network. A network PC running the Marquee software connects to the Nurse Courier system.



telecor.com

In the USA
2434 Jerauld Avenue; Niagara Falls, New York; 14305;
tel: (716) 285-8272 fax: (716) 285-8287

In Canada (Corporate)
1114 Westport Crescent, Mississauga, Ontario L5T 1G1
tel: (905) 564-0801 fax: (905) 564-0806

In the UK (International)
P.O. Box 7458,
Old Wolverton, Milton Keynes MK11 9AA Bucks, UK
tel: 44(0)1908-319228 fax: 44-1908-81-0216

All product information subject to change without notice. Copyright © 2009 by Telecor Incorporated. Telecor and the Telecor logo are Registered Trademarks of Telecor Inc.

Rev: 1.1
ETS-2702